

# PainCARE One-Touch

The ultimate in simplicity and convenience. Featuring dual channels, Conventional, Burst & Modulation modes, full digital interface, LCD screen and a comprehensive self guided user application guide.

Channels:	Dual Channel, individually isolated circuits
Amplitude:	0-80 mA
Waveform:	Asymmetrical, rectangular bi-phasic with zero DC current
PC Mode:	80 Hz, 175 $\mu$ s, 1Hr treatment time
BST Mode:	Burst of 9 pulses (200 $\mu$ s) 150 Hz. 30 Min Treatment
MOD Mode:	6 second cycle of concurrent width modulation and pulse repetition rate modulation. Width starting at 200 $\mu$ s and decreasing exponentially to 100 $\mu$ s in three seconds. Rate starting at 100 Hz, decreasing exponentially to 65 Hz and then returning to 100 Hz. 30 min treatment.
Low Battery:	Indicates if the battery goes below 6.9 volts +/- 0.2 volts.
Open Circuit:	Detects open circuit due to broken leadwires or non-contact with the skin.
Auto Off:	At the end of any program or no increase in intensity for four minutes.
Lock Out:	Button lock out with no button activity for greater than 10 seconds.
Package:	PainCARE One-Touch Unit Leadwires (2 sets) 9 Volt Battery



## AMERICAN IMEX

16520 Aston Street, Irvine CA 92606 USA  
949-553-8885 800-521-8286  
fax: 949-852-1245  
www.americanimex.com info@americanimex.com