



# PainCare

PainCare is simply efficient pain relief. Its as easy as 1,2,3.

One user-friendly device

Two Channels of 80 mA output

Three Modes of Stimulation- Conventional, Burst and Modulation

Channel: Dual, isolated between channels  
Pulse intensity: Adjustable 0-80mA peak into 500  $\Omega$  load each channel  
Pulse Rate: Pulse Rate fixed at 80Hz  
Function Modes:  
    C: Continuous mode.  
        Pulse Rate Fixed at 100Hz  
        Pulse Width Fixed at 120 $\mu$ s  
    B: Burst  
        9 pulse burst, at 100Hz, 2 bursts per second  
        Pulse Width Fixed at 250 $\mu$ s  
    M: Modulated Width.  
        Pulse Rate Fixed at 100Hz  
        Pulse Width cycles every 6 seconds from 150 $\mu$ s to 250 $\mu$ s  
  
Wave Form: Asymmetrical Bi-phasic square pulse.  
Power Source: 9-Volt Battery.  
Dimensions: 95 mm (H) x 61.5 (W) mm x 26 mm (T)  
Weight: 120 grams (battery included).



## **AMERICAN IMEX**

16520 Aston Street, Irvine CA 92606 USA  
949-553-8885 800-521-8286  
fax: 949-852-1245  
www.americanimex.com info@americanimex.com