

R-STIM II

The R-STIM II analog neuromuscular stimulator capable of contracting even the largest muscle groups. Simple analog adjustment dials for Amplitude, Rate, Width, On Ramp, Off and On times make this device a snap to set up.

| | |
|------------------|---|
| Channel: | Dual, isolated between channels |
| Pulse intensity: | Adjustable 0-100mA peak into 500 Ω load each channel, constant current |
| Wave Form: | Asymmetrical Bi-phasic square pulse |
| Pulse Rate: | 2-120 Hz |
| Pulse Width: | Fixed at 250 μ s |
| On Ramp: | Adjustable 0-8 seconds |
| Cycle ON time: | Adjustable 0-30 seconds |
| Cycle OFF time: | Adjustable 0-30 seconds |
| Cycle Modes: | Synchronous, constant or alternating |
| Power Source: | 9-Volt Battery. |
| Dimensions: | 2.31" (5.87 cm) x 3.91" (9.93 cm) x 1.00" (2.54 cm) |
| Weight: | 9.32 oz (91 grams) battery excluded |



AMERICAN IMEX

16520 Aston Street, Irvine CA 92606 USA
949-553-8885 800-521-8286
fax: 949-852-1245
www.americanimex.com info@americanimex.com